


CDC 4000

D R E A M C H A R T

Robust design, wide and bright screen for the convenient eye Vision Test environment, the widest-ranging chart and the fastest mask conversion speed are never beat virtues that makes CDC 4000 special.



Strong support for Chart functions

CDC 4000 supports 31 different kinds of general & special test charts for Vision Test along with the filtering and masking excellent capabilities.

Remote control designed by human engineering concept

Buttons are located with the full consideration of systematic & convenient operation for various types of Vision Test. One touch of a button adjusts the height of the chart to meet the customer's eye.


Flexible networking capabilities

With the implementation of high speed flexible 100K bps networking capabilities is provided in stronger and more flexible multi-to-multi communication environment supported by CAN method.

Wide & clear display coated with non-reflection method

It always provides clear image regardless of the type of lighting. It prevents excessive convergence and accomodation of the client from Vision Test.

CHART TYPE



AUTOMATIC VISUAL ACUITY TESTER - CDC 4000

| | |
|-----------------------|---|
| Refracting Distance | 5m |
| Installation Distance | 1.1 |
| Built-in Chart | 31 |
| Chart Speed | 0.45 seconds (average) |
| Masks | Single isolated, Vertical, Horizontal |
| Chart Tilting | ±20 cm |
| Lamp | 12V, 35W (Halogen) |
| Automatic Power Off | 10 minutes |
| Dimension | 400mm(W)x280mm(D)x545mm(H) |
| Weight | 20kg |
| Power Supply | AC 100~120/200~240, 50 / 60Hz, 2.0/1.0A |
| Power Consumption | 90W |

Designs and details can be changed without prior notice for the purposes of improvement

Huvitz

Huvitz Building 689-3, Geumjeong-dong
Gunpo-si, Gyeonggi-do, 435-862, Korea

Tel : +82-31-442-8868
Fax : +82-31-477-8617
[http : //www.huvitz.com](http://www.huvitz.com)

Distributed by