

Huvitz Instrument Table



CIT 4000 Instrument Table

Huvitz Instrument Table



Free load for Huvitz Instruments

You can easily install Huvitz Instruments, Slit Lamp HS-5000 or Auto Refractometer.


Easy up, down movement by simple pedal operation

INSTRUMENT TABLE - CIT 4000


SPECIFICATIONS

Table Top for Slit Lamp HS-5000


(Table Top can be replaced as below to install Huvitz Slit Lamp HS-5000)



Front View



Side View



Slit Lamp Table Top View

Dimensions	Same as above pictures
Weight	26kg(including Table Top for Auto Refracmetor)
	29kg(including Table Top for Slit Lamp HS-5000)
Power Supply	AC 220-230V, 50/60Hz
Travel Length	250 mm
Speed	22mm/sec(under 40kg), 11mm/sec(under 60kg)
Load	Max 60kg
Others	Up, down movement operation pedal included

Designs and details can be changed without prior notice for its improvements.



Huvitz Building 689-3, Geumjeong-dong

Tel : ++82-31-442-8868

Fax : ++82-31-477-8617

http : //www.huvitz.com

Distributed