




Seeing is **believing.**



VOLK **iNview**

iPhone fundus camera



Get with the **iN** crowd.

Leverage the power and convenience of the Apple iPhone with the trusted quality of Volk optics.

Quickly & effortlessly capture fundus images for visualization & patient education.

Helps facilitate patient discussions related to disease progression and treatment plan.

- > Free mobile application available in the Apple App Store (search Volk iNview)
- > Offers 1 Megapixel resolution with a static 50° field of view
- > View the peripheral retina dynamically out to 80°
- > Available auto-capture & forced-capture imaging modes
- > Mydriatic; requires minimum 5mm pupil
- > HIPAA-compliant storage and export from iPhone to PC or Mac
- > Compatible with Apple iPhone 6s/6/5s and iPod Touch (Gen 6)