



The future of vision

ANICAN



ANICAN Capsular Tension Ring

Model	ACV				
Reference	ACV001002	ACV001102	ACV001202	ACV001302	ACV001402
Open diameter	10 mm	11 mm	12 mm	13 mm	14 mm
Closed diameter	9 mm	10 mm	10 mm	11 mm	12 mm
Brand name	ANICAN CAPSULAR TENSION RING				
Product description	ANICAN Capsular Tension Ring is an ophthalmic device designated to stabilize the capsule in canine cristalline extraction procedures.				
Indications	The purpose of these implants is to expand the capsule in which the crystalline is located so that an intraocular lens can be inserted.				
Material	Medical grade Polymethyl Methacrylate (PMMA) with UV filter.				
Design	AC: round section.				
Sterilization method	Ethylene Oxide (ETO).				
Supply	Single and sterile in an easy-handling carrier.				
Shelf-life	4.5 years				

REV: 01