




The future of vision

BLUEVET



BLUEVET

Trypan Blue Ophthalmic Solution

Reference	BLUEVET
Model	Trypan blue ophthalmic solution at 0.06% W/V.
Brand name	BLUEVET ophthalmic solution.
Product description	BLUEVET is an intraocular and biocompatible staining solution for use in anterior segment surgery. The use of BLUEVET provides a good visual image of the capsular bag. By staining the cristalline's anterior chamber, its visibility is improved, which means the capsulorhexis is simplified and the risk of tearing is minimized.
Composition	1ml of Bluevet contains 0.6mg of Trypan Blue Ophtahmic Solution.
Indications	BLUEVET is indicated in eyes with mature cataracts, poor fundus reflex or small pupils.
Sterilization method	Steam.
Supply	Box: 5 units. Contains a 1ml prefilled glass syringe and Luer lock 27G cannula.
Shelf-life	3 years

REV: 01