



The future of vision

SIOBALVET



SIOBALVET Silicone Oil

| | |
|-----------------------------|--|
| <i>Reference</i> | S5000VET |
| <i>Brand name</i> | SIOBALVET |
| <i>Product description</i> | Tamponade agent for retinal detachment. |
| <i>Material</i> | High purity silicone oil free of low molecular weight components. |
| <i>Viscosity</i> | 5 000 cSt |
| <i>Indications</i> | <ul style="list-style-type: none">• Tamponade agent in complex cases of retinal detachment.• Silicone oil is used as a vitreous alternative.• For veterinary use only. |
| <i>Sterilization method</i> | Steam. |
| <i>Supply</i> | Single and sterile in a 10ml prefilled glass syringe. |
| <i>Shelf-life</i> | 3 years |

REV: 01