

# REFRACTION UNIT TRU-800

SMART REFRACTION UNIT

## DELIGHT IN SIGHT

Smart solution. For everybody.



- Easy to install (right or left)
- Low shipment costs
- Space saving (1.08 m<sup>2</sup>)
- Incl. chair SC-1000
- Magnetic table lock
- Timeless and modern design

 **TOMEY**  
TECHNOLOGY AND VISION

# THE TOMEY REFRACTION UNIT TRU- 800



## QUALITY IN DETAIL

Perfect if you need a tiny and clever solution for your refraction environment. Much smaller than standard refraction units does it even fit to very small rooms. For a proper refraction just attach all essential instruments as manual/automatic phoropter and chart projector to the TRU-800. If you want space for hand-held units or trial lenses - the optional trial lens cabinet TC-1000 meets with the design.



Rotatable control board panel



Phoropter arm PA-2000 and PA-1000 (optional)



Trial lens cabinet TC-1000 (optional)

## SPECIFICATIONS

**Seat height** 42 cm (SC-1000)

**Seat elevation** 42 - 57 cm

**Nominal lifting capacity** max. 120 kg

**Projector column height** 175 cm

**Overhead light** white LED 3.5 V

**Table plate** 44 x 87 cm

**Table lock** magnetic

**Table height** 81 cm

**Weight** approx. 150 kg


**Power consumption** 300 VA

**Power supply** 230 V

**Connectors** 2 x power outlet (230V),  
1 x 5-pin (6V/12V)

### Optionals:

Phoropter arm, trial lens cabinet 61 x 61 x 79 cm (WxDxH),  
ophthalmoscope holder



Dimensions

2016/08 - subject to change without notice

**TOMEY EUROPE**  
TOMEY GmbH  
Am Weichselgarten 19a  
91058 Erlangen, Germany

Phone +49 9131 777 10, Fax +49 9131 777 1 20  
Email info@tomey.de

**TOMEY ASIA-PACIFIC**  
TOMEY CORPORATION JAPAN  
2-11-33 Noritakeshinmachi  
Nishi-ku, Nagoya 451-0051, Japan

Phone +81 52 581 5327, Fax +81 52 561 4735  
Email intl@tomey.co.jp

**TOMEY**  
TECHNOLOGY AND VISION

www.tomey.de