



The future of vision

AIALA DRY



AIALA DRY

Intraocular Lenses

Brand name	AIALA DRY	AIALA DRY PRELOADED
References	LLASHP60	LLASHP60-PL
Description	Aspherical hydrophobic intraocular lens. One piece design. Posterior chamber.	
Design	C-LOOP, antistick haptics. Square Edge (360°)	
Material	Hydrophobic acrylic, GF, with filter for UV and filter for blue light.	
Optical diameter	6.0 mm	
Total length	13.00 mm	
Angle	0°	
Diopters range	-5 D to +40 D in steps of 1 D / +10 D to +27 D in steps of 0.5 D	
Water content	<4%	
Constant A / ACD	A-constant (optical): 119.7 A-constant (U/S): 119.2 ACD: 5.66	
Sterilization method	Ethylene Oxide (ETO).	
Supply	Single unit, sterile.	
Shelf-life	3 years	

* Product with CE marking

* Pre-loaded intraocular lens injector
Ref.: LLASHP60-PL

REV: 02. 06/2019