




The future of vision

AJL BBG



# AJL BBG Ophthalmic Solution

Brand name	AJL BBG
Reference	AJL BBG 0.025%
Description and Indications	<p>AJL BBG 0.025% is an ophthalmic sterile, non-pyrogenic, isotonic solution, ready for use in vitreoretinal surgery. The BBG solution is capable of effectively staining the inner limiting membrane (ILM).</p> <p><u>Indications:</u></p> <p>Greater affinity with the inner limiting membrane (ILM) and lesser with the epiretinal membrane, allowing an improved view of the ILM and facilitating the peeling in surgical operations for macular diseases, such as Macular Hole or Macular Fibrosis.</p> <p><u>Properties:</u></p> <ul style="list-style-type: none"> <li>• Ready-to-use sterile solution: does not require mixing or filtering.</li> <li>• Excellent biocompatibility.</li> <li>• No adverse side-effects.</li> <li>• Easy elimination of residual product.</li> </ul>
Composition	Brilliant blue G sterile ophthalmic solution at 0.025%
Sterilization method	Steam.
Supply	Sterile, single unit. Contains a 1 ml prefilled glass syringe.
Shelf-life	3 years

\* Product with CE marking

REV: 04. 06/2019