



ESNOOPER V2000



The future of vision

ESNOOPER Glaucoma Implant

Brand name	ESNOOPER
Reference	V2000
Description e indication	The ESNOOPER V2000 scleral implant is a non-absorbable implant that provides a fixed, permanent intrascleral space that facilitates the drainage of aqueous humor. It is indicated for non-penetrating deep sclerectomy (NPDS) in open-angle glaucoma surgery. The ESNOOPER is placed on the scleral bed, this method favors the trabecular way, which is the classic drainage way of the aqueous humor
Dimensions	3.00mm X 1.40mm X 2.85mm.
Orifices	2 (suture/drainage).
Longitudinal orifice	Drainage.
Material	Hydroxyethyl Methacrylate (HEMA).
Sterilisation method	Steam.
Package	The ESNOOPER V2000 implant is supplied single and sterile in a vial, housed by a clamp support that is immersed in high purified water to maintain it properly hydrated.
Shelf-life	3 years
* Product with CE marking	

NPDS surgical instruments (recommended)

Brand name	KIT EPNP	Product code	AJL- KIT-001
Definition	NPDS (non-penetrating deep sclerectomy) Kit .		
Description	<p>The NPDS Kit is presented non sterile in a rigid plastic blister with the following items (*):</p> <ul style="list-style-type: none"> • EPNP-2001 Liebermann Speculum 14mm. • EPNP-2002 Double Scleral Marker. • EPNP-2003 Mermoud Spatula. • EPNP-2004 Mermoud Clamp. • EPNP-2005 Curved micro-scissor. • PST-4600A Sterilisation box (10x16x2 cm). 		
* The instruments can be supplied as a set (NPDS Kit) or separately.			



REV: 03. 04/2019