



KSL-H  
Keeler

---

Top down illumination



## KSL-H3

### Analog Slit Lamp

Tower Illumination

3 Step Drum  
(10x, 16x, 25x)

## KSL-H5

### Analog Slit Lamp

Tower Illumination

5 Step Drum  
(6x, 10x, 16x, 25x, 40x)

## KSL-H3-DR

### Digital Ready Slit Lamp

Available Add-Ons:

- TTI Beam Splitter
- Canon SLR Camera
- AOS Slit Diagnostic Software
- Kapture Software

Tower Illumination

3 Step Drum  
(10x, 16x, 25x)

## KSL-H5-DR

### Digital Ready Slit Lamp

Available Add-Ons:

- TTI Beam Splitter
- Canon SLR Camera
- AOS Slit Diagnostic Software
- Kapture Software

Tower Illumination

5 Step Drum  
(6x, 10x, 16x, 25x, 40x)

## KSL-H3-D

### Digital Slit Lamp

2 Megapixel:  
1600 x 1200 Camera

Kapture Software

Available Add-On:

- AOS Slit Diagnostic Software

Tower Illumination

3 Step Drum  
(10x, 16x, 25x)

## KSL-H5-D

### Digital Slit Lamp

2 Megapixel:  
1600 x 1200 Camera

Kapture Software

Available Add-On:

- AOS Slit Diagnostic Software

Tower Illumination

5 Step Drum  
(6x, 10x, 16x, 25x, 40x)