


DESIGNED WITH CONVENIENCE IN MIND

Ahmed® ClearPath™ - the newest valveless option offered from the trusted Ahmed brand of glaucoma drainage devices.

Model 350

Plate surface is designed to avoid muscle attachment points




Contour designed for fit between the implant and globe

A pre-threaded 4-0 polypropylene ripcord is included with each model

Suture points are positioned anteriorly providing easier access for implantation

Model 250

True single-quadrant implanted between the muscles



23 gauge needle included with every model for added convenience



Specifications

	Ahmed® ClearPath Model 250	Ahmed® ClearPath Model 350	Ahmed® Glaucoma Valve Model FP7	Ahmed® Glaucoma Valve Model FP8
Type	Valveless	Valveless	Valve	Valve
Surface Area	250 mm ²	350 mm ²	184 mm ²	102 mm ²
Plate Length	16.7 mm	16.5 mm	17.0 mm	11.0 mm
Plate Width	18.5 mm	30.5 mm	13.0 mm	11.0 mm
Plate Thickness	0.8 mm	0.8 mm	2.1 mm	2.1 mm
Tube Length	32 mm	32 mm	25.4 mm	25.4 mm
Fenestrations	4	4	3	0
Material	Barium-impregnated silicone		Medical-grade silicone	

Reimbursement

CPT® Code	Description of Procedure
66179	Aqueous shunt to extraocular reservoir; without patch graft
66180	Aqueous shunt to extraocular reservoir; with patch graft

For more information, contact your local representative.



New World Medical
10763 Edison Court
Rancho Cucamonga, CA 91730
(909) 466-4304