



The future of vision

LIOFELINE



LIOFELINE

Feline Intraocular Lenses

Brand name	LIOFELINE		
Model	CAT70		
Reference	CAT701200	CAT701300	CAT701400
Optic diameter	7.0 mm	7.0 mm	7.0 mm
Overall diameter	12.00 mm	13.00 mm	14.00 mm
Diopters	53.5 D	53.5 D	53.5 D
Descripción	Liofeline Intraocular Lens is an ophthalmic device designed for implantation in canine crystalline extraction procedures. Liofeline intraocular lens has been specifically designed for most canine breeds, establishing different lens diameters for a better adaptation of the lens within the capsular bag.		
Material	Hydrophilic Acrylic with UV filter and an additional filter that blocks blue light, thus offering additional protection to the retina.		
Optic design	Biconvex / Aspherical	Biconvex / Aspherical	Biconvex / Aspherical
Sterilization method	Steam.		
Supply	Sterile. Single unit. Implantation with Medice Accuject™ 3.0mm - 1P injector set is recommended.		
Shelf-life	3 years		

REV: 01. 03/2021