

Precision Ophthalmic Instruments Single Use

ACROfine[®] dispotz



www.eyeline.org



About Eyeline

Eyeline Organisation Ltd takes pride in presenting single-use instruments under the trademarks Acrofine® dispøz® and Derafine dispøz® for our packs, specifically designed for minor operations in primary care and dermatologists. Designed and manufactured to the highest standards by surgeons for surgeons, these instruments are made from the same high-grade surgical stainless steel as reusable instruments, and then hand finished to perfection.

They feel and perform exactly like their reusable counterparts.

Instruments are checked by eye surgeons themselves and each and every instrument is quality controlled to the highest standard.

At Eyeline we believe firmly that single-use instruments are the way forward, not only for ophthalmic surgery but for many other surgical sub-specialties, so we are continually developing new procedure packs to make life easier and more reliable for both surgeons and hospitals. We provide a wide range, however some instruments may be available by special order only.

Ophthalmic instruments need careful handling. Expensive reusable instruments must be cared for every time they are processed. Sterilisation at super-centres takes away some of the control hospitals have over this care. Mass processing increases the risk of damage.

Defects in reprocessed instruments may only become obvious at the critical moment during surgery. Single-use procedure packs guarantee a perfect set of instruments every time.

PROCEDURE PACKS

Acrofine® dispøz® & Derafine dispøz® procedure packs are Eyeline's flagship products. Based on the concept of a tray of instruments, they offer all the essentials for the surgical procedure in a simple-to-use, transparent pack.



Phaco cataract procedure pack 9990-6D

Contents

9960D-Lid speculum
9989-PD-Swab holder x2
9929D-Conjunctival scissors straight pointed
9937D-Needle bender
9906D-Corneal forceps toothed with platform
9910D-Ultrata curved forceps
9945D-Sinskey dialler
9977D-Cannula 21G
9981D-Hydrodissection cannula (Helsinki style)



Phaco Lite procedure pack 9804D

Contents

9960D-Lid speculum
9910D-Ultrata curved forceps
9906D-Corneal forceps toothed with platform
9948D-Mushroom micro manipulator (Barret style)

Variants (Mushroom micro manipulator replaced)

9804-1D with Nagahara chopper
9804-2D with Micro spatula (Tayyib)
9804-3D with Green chopper RH surgeon
9804-4D with Micro finger
9804-5D with Micro paddle (Drysdale)



Phaco extension procedure pack 9805D

Contents

9908D-McPherson forceps
9945D-Lens dialler Sinskey style

Variants (Sinskey lens dialler replaced)

9805-1D with Nagahara chopper
9805-2D with Micro spatula (Tayyib)
9805-3D with Green chopper RH surgeon
9805-4D with Micro finger
9805-5D with Micro paddle (Drysdale)
9805-7D with Mushroom micro manipulator (Barret)

PROCEDURE PACKS

Just open one tray for one operation and it's all there. Acrofine® & Dermafine dispaz® instruments and procedure packs are the way forward both for the modern operating theatre and the community-based surgeon. Eyeline's procedure packs consist of:



Intravitreal Lite procedure pack 9995D
Variants: 9995-1D with lid speculum (phaco) instead of Acrofine guarded speculum




Intravitreal injection procedure pack 9993D




Sub-Tenons LA Lite procedure pack 9999D



Sub-Tenons LA admin. procedure pack 9991D



Chalazion procedure pack 9992D



Suture removal packs 9800D(left) & 9801D(right)



Micro suturing procedure pack 9996D



Syringing & probing procedure pack 9994D



Neonatal examination pack 9998D



Lid lesion biopsy pack 9803D



Dermafine procedure pack 9807D



Oculoplastic procedure pack 9809D

FORCEPS

Acrofine® dispøz® forceps are available in a wide range, designed for all branches of ophthalmic surgery from oculoplastics to intraocular work. The tips are hand finished to ensure the highest quality.



Corneal
plain
forceps
9900D



Corneal
tying
forceps
platform length
6mm
9902D



Corneal
notched
forceps
with platform
notch
diameter
0.20mm
9905D



Corneal
toothed
forceps with
platform
tooth height
0.12mm
9906D



McPherson
forceps
angle to tip
8mm
9908D



Ultracurved
forceps
angle to tip
13mm
9910D



St. Martins
forceps
tooth
length
0.27mm
9917D



St. Martins
forceps
large
9909D



Jayles
forceps
tooth
length
0.27mm
9918D



Moorfields
forceps
serrated tip
9919D



Adson's
plain
forceps
E1000D



Adson's
toothed
forceps
E1002D

MICRO MANIPULATORS

We understand the critical importance of 'second instruments' in phaco surgery. Acrofine® dispaz® micro manipulators are again hand finished. A wide range of tips are available to cater for the preferences of different surgeons.



Micro
paddle
(Drysdale)
9939D



Micro
spatula
(Tayyib)
tip length
1mm
9940D



Micro
finger
tip length
2mm
9941D



Chopper
Nagahara
style
tip length
1.25mm
9942D



Chopper
Green style
right handed
surgeon
tip length
2mm
9943D



Chopper
Green style
left handed
surgeon**
tip length
2mm
9944D



Lens
dialler
Sinskey
style
tip length
0.9mm
9945D



T-dialler
Lester
style
tip length
1.4mm
9946D



Y-dialler
Bechert
style
tip length
0.5mm
9947D



Mushroom
Barret
style
tip thickness
0.65mm
9948D



Iris
repositor
flat blade
9949D



Iris
repositor
round
blade
9959D

SCISSORS / NEEDLE HOLDERS / SPECULA

Acrofine® dispøz® ophthalmic scissors have been carefully selected for a whole range of uses from cutting bandages to managing very fine intraocular tasks. Being single-use instruments they are always sharp. The Acrofine® dispøz® needle holder has been designed to combine the finesse of a micro needle holder and fine suture tier with the utility and feel of an instrument capable of handling larger needles used in lid surgery.



Vannas
curved
scissors
blade length
8mm
9920D



Vannas
straight
scissors
blade length
8mm
9921D



Williamson
Noble
scissors
cured pointed
blade length
18.2mm
9924D



Williamson
Noble
scissors
curved blunt
blade length
18.2mm
9925D



Westcott
scissors
blunt
9927D



Conjunctival
curved
blunt
scissors (bow)
cut length
24mm
9928D



Conjunctival
straight pointed
scissors (bow)
cut length
24mm
9929D



Enucleation
scissors
cut length
26mm
9931D



Dressing
scissors
cut length
40mm
9932D



Stevens
scissors
(supercut
blades)
E2000D



Needle
holder
(Bates)
jaw length
10mm
9935D



Needle
holder
(Castroviejo
with lock)
9936D



Halsey
Webster
needle
holder
E3000D



Lid speculum
(adult)*
standard
9960D
& temporal
access
9961D



Lid
speculum -
Acrofine
guarded
9964D



Wire lid
retractor
(Cox Jaffe)
9963D



Desmarres
retractor
12mm
E4010D



Cats paw
retractor
(Blair-Rollet)
8mm
E4020D



Squint
hook
(adult-Graefe)
E4030D



Skin hook
(Gilles)
E4000D

CANNULAE / LID / MISCELLANEOUS



Braunstein
calliper
3.5/4.00mm
9965D



Vectis /
scleral
indenter
9966D



Lid clamp
(chalazion)
ring
diameter
10.5mm
9970D



Double-ended
curette
(chalazion)
spoon size
1.7mm &
2.6mm
9972D



Scalpel
handle no.3
with no.11 blade
9973-PD



Barron's
scalpel
handle with
no.15 blade
E5000-PD



Punctum
dilator
(Nettleships)
9974D



Punctum
dilator /
seeker
(double-
ended)
9976D



Lacrimal
probes
(set of 3)
diameters
0.5, 0.7
& 1.0mm
9975D



Irrigation /
aspiration
cannula
Simcoe
style
port 0.4 mm
9985D



Sub-
Tenons
cannula*
9979D



Hydro-
dissection
cannula
angled*
25G
9980D



Hydro-
dissection
cannula
Helsinki
style*
26G
9981D



Hydro-
dissection
cannula
Bolger style*
26G
9982D



Towel
clip
9988D



Artery
forceps
curved
9986D



Artery
forceps
straight
9987D



Swab
holder
forceps
9989-PD

Product Table

Instrument	Cat. No.	Instrument	Cat. No.	Instrument	Cat. No.
PROCEDURE PACKS		FORCEPS		Sinsky style Implant lens dialler	9945D
Phaco cataract procedure packs:		Corneal plain forceps	9900D	T-dialler Lester style	9946D
with Sinsky lens dialler	9990-6D	Corneal tying forceps	9902D	Y-dialler Bechert style	9947D
Phaco cataract Lite packs:		Corneal notched forceps with platform	9905D	Mushroom micromanipulator Barrett style	9948D
with mushroom micro manipulator	9804D	Corneal toothed forceps with platform	9906D	Iris repositor-flat blade	9949D
with Nagahara chopper	9804-1D	McPherson forceps	9908D	Iris repositor-round blade	9959D
with micro spatula (Tayyib)	9804-2D	Utrata curved forceps	9910D	SPECULA	
with Green chopper RH surgeon	9804-3D	St. Martins forceps	9917D	Lid speculum (phaco)	9960D
with micro finger	9804-4D	St. Martins forceps large	9909D	Lid speculum (temporal access)	9961D
with micro paddle (Drysdale)	9804-5D	Jayles forceps	9918D	Lid speculum (neonatal)	9962D
Phaco extension procedure pack		Moorfields forceps	9919D	Wire lid retractor, Cox Jaffe (packed in pairs)	9963D
with Sinsky dialler	9805D	Adson's plain forceps	E1000D	Lid speculum - Acrofine guarded	9964D
with Nagahara chopper	9805-1D	Adson's toothed forceps	E1002D	MISCELLANEOUS	
with Micro spatula (Tayyib)	9805-2D	SCISSORS		Braunstein calliper (plastic)	9965-PD
with Green chopper RH surgeon	9805-3D	Vannas curved scissors	9920D	Vectis / scleral indenter	9966D
with Micro finger	9805-4D	Vannas straight scissors	9921D	Lid clamp (chalazion)	9970D
with Micro paddle (Drysdale)	9805-5D	Williamson Noble curved pointed scissors	9924D	Double ended currette (chalazion)	9972D
with Mushroom micro manipulator (Barrett)	9805-7D	Williamson Noble curved blunt scissors	9925D	Scalpel handle no.3 with no.11 blade & guard	9973-PD
Sub-Tenons LA admin. procedure pack	9991D	Westcott scissors blunt	9927D	Barron's scalpel handle with no.15 blade & guard	E5000-PD
Sub-Tenons Lite LA admin. procedure pack	9999D	Conjunctival curved blunt scissors (bow)	9928D	Punctum dilator (Nettleship's)	9974D
Chalazion procedure pack	9992D	Conjunctival straight pointed scissors (bow)	9929D	Punctum dilator / seeker (double-ended)	9976D
Intravitreal injection procedure pack	9993D	Enucleation scissors	9931D	Lacrimal probes (set of 3)	9975D
Intravitreal injection Lite procedure pack (guarded speculum)	9995D	Dressing scissors	9932D	Sub-Tenons cannula*	9979D
Intravitreal injection Lite procedure pack (phaco speculum)	9995-1D	Stevens scissors (supercurt blades)	E2000D	Hydrodissection cannula (angled)*	9980D
Syringing & probing procedure pack	9994D	NEEDLE HOLDERS		Hydrodissection cannula (Helsinki style)*	9981D
Micro suturing procedure pack	9996D	Needle holder (Bates)	9935D	Hydrodissection cannula (Bolger style)*	9982D
Neonatal examination procedure pack	9998D	Needle holder (Castroviejo with lock)	9936D	Irrigation / aspiration cannula (Simcoe)	9985D
Suture removal pack (Vannas scissors)	9800D	Needle holder (Halsey Webster)	E3000D	Artery forceps curved 117 mm	9986D
Suture removal pack (Williamson Noble scissors)	9801D	MICRO MANIPULATORS		Artery forceps straight 117 mm	9987D
Lid lesion biopsy pack	9803D	Micro paddle (Drysdale)	9939D	Towel clip	9988D
Minor operations / Dermatological pack	9807D	Micro spatula (Tayyib)	9940D	Swab holder (packed in pairs)	9989-PD
Oculoplastic procedure pack	9809D	Micro finger	9941D	Skin hook (Gilles)	E4000D
		Chopper Nagahara style	9942D	Desmarres retractor	E4010D
		Chopper Green style (right handed surgeon)	9943D	Cats paw retractor (Blair-Rollet)	E4020D
		Chopper Green style (left handed surgeon)**	9944D	Squint hook (adult-Graefe)	E4030D

* Not available separately

** Special order only

www.eyeline.org

