



## All Pupil II

The All Pupil II is based on the design of the number one selling All Pupil indirect. Look out for the new features designed to make your examinations easier and more comfortable.

**Go Wireless patented technology**

Leaves you free to move around without any wire or cord constraints.

**Outstanding brightness for difficult examinations**

With 2000 lux of light, which can be turned down to 2% of total illumination when necessary, the All Pupil II offers maximum illumination with total control.

**Lightweight and comfortable**

Weighing just 520 grams with extra padded materials.

**Simple controls**

The All Pupil II has a unique single control to change aperture size and position in one movement to view through all pupil sizes.

**IR and UV barrier**

There is a built in full time IR and UV barrier which guarantees improved patient protection and safety.

**New All Pupil II LED light source**

Plugs simply into the back of the All Pupil II for brighter, whiter illumination and freedom from purchasing bulbs. The LED life expectancy is up to 10,000 hours and it has a 5 year warranty.



**Part numbers**

1204-P-3041	Wireless All Pupil II Kit A with Power / Accessories
1204-P-3043	Wireless All Pupil II with Lithium Battery and Holder
1204-P-3038	Wired All Pupil II Kit A with Power / Accessories
1204-P-3036	Wired All Pupil II without Power

**Accessories**

1205-P-5014	All Pupil II LED
1202-P-7117	All Pupil II Teaching Mirror
1012-P-7003	All Pupil II/All Pupil/Vantage Bulbs (pack of 2)
3412-P-7000	All Pupil II Carrying Case
1201-P-6075	Small Depressor
1201-P-6067	Large Depressor

For Wired Power Supplies see page 13

