

Ophthalmic Surgical Chair

SC-1

Ergonomics, functionality, comfort



SC-1 surgical table has been designed to ensure maximum comfort for the patients and surgeons during ophthalmic surgical procedures with complete seat tilt, Trendelenburg position, and different accessories available.

Full range of patient positioning from sitting to lying position gives unprecedented opportunity to position patients at any age in a comfortable way.

A unique design of the patient's head rest gives the surgeon precise control over the patient's head during surgical procedures.

A robust construction of the base with special casters allow the surgical table to be freely moved in and out the operation theatre with the patient sitting on.

Central locking gear assures stability during surgery.

Battery/mains operated controls assure maximum flexibility of the system.

Wide selection of upholstery colours available.

Technical specification:

Overall length:	1750mm-1850mm (flat position)
Overall width:	925mm
Net weight:	100kg
Max. load of table:	250kg

Adjustment range:

Height	610 - 750mm
Trendelenburg position	15°
Detachable arm supports	✓
Detachable cushioning	✓

Accessories:

Anaesthesia support/goose neck tube holder	✓
Push-handles	✓

Electrical data:

Operating voltage	24 VDC
Battery capacity	12Ah
Operating cycle at full load	40 cycles

Available colours:

