



AJL
ophthalmic

ALOS



Implantes Palpebrales

Nombre de marca	ALOS
Modelos	ORO (TH - G) y PLATINO (PTH - P)
Descripción e indicaciones	• Materiales biocompatibles. • Debido a su diseño esférico y fino, los implantes se adaptan a la curvatura del globo ocular para conseguir un mejor resultado estético. • Bordes redondeados para evitar la extrusión del implante a largo plazo. • Orificios de sutura con canales de fijación. • Pureza del 99,9%(Platino) y 99,9%(Oro). • Estéril. <u>Indicaciones:</u> Los implantes palpebrales están diseñados para la corrección del lagofthalmos.



Nombre de marca	ALOS SIZER SET
Referencias	SET MEDIDOR PESAS TH SET MEDIDOR PESAS G
Peso (g)	0.6-0.8-1.0-1.2-1.4-1.6

ORO																		
Referencias	TH 06	TH 08	TH 10	TH 12	TH 14	TH 16	TH 18	G 06	G 08	G 10	G 12	G 14	G 16	G 18	G 20	G 22	G 24	G 26
Peso (g)	0.6	0.8	1.0	1.2	1.4	1.6	1.8	0.6	0.8	1.0	1.2	1.4	1.6	1.8	2.0	2.2	2.4	2.6
Longitud (mm)	10.5	12.5	15.5	17.5	19.5	21.5	24.5	7.8	10.3	12.3	13.9	15.4	16.9	18.9	20.5	20.8	22.3	23
Ancho (mm)	5.7	5.7	5.7	5.7	5.7	5.7	5.7	5	5	5	5	5	5	5	5	5	5	5
Grosor (mm)	0.6	0.6	0.6	0.6	0.6	0.6	0.6	1.1	1.1	1.1	1.1	1.1	1.1	1.1	1.1	1.1	1.1	1.1
Método de esterilización	Vapor																	
Caducidad	3 años																	
PLATINO																		
Referencias	PTH06	PTH08	PTH 10	PTH 12	PTH 14	PTH 16	PTH 18	P 06	P 08	P 10	P 12	P 14	P 16	P 18				
Peso (g)	0.6	0.8	1.0	1.2	1.4	1.6	1.8	0.6	0.8	1.0	1.2	1.4	1.6	1.8				
Longitud (mm)	10.2	12.4	15.0	17.0	18.8	21.0	24.0	7.6	10	11.9	13.5	15.0	16.3	17.4				
Ancho (mm)	5.7	5.7	5.7	5.7	5.7	5.7	5.7	5	5	5	5	5	5	5				
Grosor (mm)	0.6	0.6	0.6	0.6	0.6	0.6	0.6	1.1	1.1	1.1	1.1	1.1	1.1	1.1				
Método de esterilización	Vapor																	
Caducidad	3 años																	

* Producto con marcado CE



REV: 05. 07/2022



ALOS



Palpebral Implants

Brand name	ALOS
Models	GOLD (TH - G) and PLATINUM (PTH - P)
Description and Indications	• Biocompatible materials • Due to their fine, spherical design, the implants adapt to the curvature of the eyeball to achieve a better aesthetic result. • Rounded edges to avoid the extrusion of the implant in the long-term. • Suture holes with fixation channels. • Purity of 99,9%(Platinum) and 99,9%(Gold). • Sterile. <u>Indication:</u> Eyelid implants are designed for the correction of lagophthalmos.



Brand name	ALOS SIZER SET
References	SET MEDIDOR PESAS TH SET MEDIDOR PESAS G
Weight (g)	0.6-0.8-1.0-1.2-1.4-1.6

GOLD

References	TH 06	TH 08	TH 10	TH 12	TH 14	TH 16	TH 18	G 06	G 08	G 10	G 12	G 14	G 16	G 18	G 20	G 22	G 24	G 26
Weight (g)	0.6	0.8	1.0	1.2	1.4	1.6	1.8	0.6	0.8	1.0	1.2	1.4	1.6	1.8	2.0	2.2	2.4	2.6
Length (mm)	10.5	12.5	15.5	17.5	19.5	21.5	24.5	7.8	10.3	12.3	13.9	15.4	16.9	18.9	20.5	20.8	22.3	23
Width (mm)	5.7	5.7	5.7	5.7	5.7	5.7	5.7	5	5	5	5	5	5	5	5	5	5	5
Thickness (mm)	0.6	0.6	0.6	0.6	0.6	0.6	0.6	1.1	1.1	1.1	1.1	1.1	1.1	1.1	1.1	1.1	1.1	1.1
Sterilization method	Steam																	
Shelf-life	3 years																	

PLATINUM

References	PTH06	PTH08	PTH 10	PTH 12	PTH 14	PTH 16	PTH 18	P 06	P 08	P 10	P 12	P 14	P 16	P 18
Weight (g)	0.6	0.8	1.0	1.2	1.4	1.6	1.8	0.6	0.8	1.0	1.2	1.4	1.6	1.8
Length (mm)	10.2	12.4	15.0	17.0	18.8	21.0	24.0	7.6	10	11.9	13.5	15.0	16.3	17.4
Width (mm)	5.7	5.7	5.7	5.7	5.7	5.7	5.7	5	5	5	5	5	5	5
Thickness (mm)	0.6	0.6	0.6	0.6	0.6	0.6	0.6	1.1	1.1	1.1	1.1	1.1	1.1	1.1
Sterilization method	Steam													
Shelf-life	3 years													

* Product with CE marking



REV: 05. 07/2022