




AIALA

IOLs

Referencia	F601250
Nombre de marca	AIALA
Descripción	Lente intraocular hidrofílica. Diseño monobloque. Cámara posterior.
Diseño	Hápticas C.
Material	Acrílico hidrofílico
Filtro	UV-A
Diámetro de la óptica	6.0 mm.
Longitud total	12.50 mm.
Angulación	0°
Rango de dioptrías	-6 D a +40 D en incrementos de 1 D / +10 D a +25 D en incrementos de 0,5 D
Borde	Cuadrado (360°)
Contenido de agua	25%
Constante A / ACD	Constante A: (óptica) 118 Constante A: (u/s) 117.3 ACD: 4.9
Método de esterilización	Vapor.
Presentación	Envase unitario, estéril.
Sistema inyector recomendado	Medicel.
Caducidad	3 años.

* Producto con marcado CE



REV: 05. 07/2022



AIALA



IOLs

Reference	F601250
Brand name	AIALA
Description	Hydrophilic intraocular lens. One piece design. Posterior chamber.
Desing	C Haptics.
Material	Hydrophilic acrylic.
Filter	UV-A
Optical diameter	6.0 mm
Total length	12.50 mm
Angle	0°
Diopters range	-6 D to +40 D in steps of 1 D / +10 D to +25 D in steps of 0.5 D
Edge	Square (360°)
Water content	25%
Constant A / ACD	Constant A: (optical) 118 Constant A: (u/s) 117.3 ACD: 4.9
Sterilization method	Ethylene Oxide (ETO).
Supply	Single unit, sterile.
Injector system recommended	Medicel.
Shelf-life	3 years.

* Product with CE marking



REV: 05. 07/2022