



PROTECTOR OCULAR

Protector Ocular

Nombre de marca	AJL PROTECTOR OCULAR		
Referencias	PROCUL-S-AJL	PROCUL-M-AJL	PROCUL-L-AJL
Tamaño	Pequeño: 22 x 19 x 1.7 mm	Mediano: 24 x 21 x 1.7 mm	Grande: 25 x 22 x 1.7 mm
Descripción e indicaciones	El PROTECTOR OCULAR AJL es un dispositivo médico especialmente diseñado para proteger el globo ocular. Está diseñado para proteger el globo ocular durante la cirugía de párpados.		
Material	Acrílico de alta resistencia y acabado superficial apto para el uso ocular.		
Método de esterilización	Óxido de etileno (ETO).		
Presentación	Envase unitario, estéril. Se incluye una ventosa de succión, sin látex, que facilita la colocación y extracción del protector ocular.		
Caducidad	4,5 años		

* Producto con marcado CE



REV: 05. 07/2022



OCULAR SHIELD

Ocular Shield

Brand name	AJL OCULAR SHIELD		
References	PROCUL-S-AJL	PROCUL-M-AJL	PROCUL-L-AJL
Size	Small: 22 x 19 x 1.7 mm	Medium: 24 x 21 x 1.7 mm	Large: 25 x 22 x 1.7 mm
Description and Indications	AJL OCULAR SHIELD is a medical device specially designed to protect the eyeball. It's designed to protect the eyeball during eyelid surgery.		
Material	High resistance acrylic material suitable for ocular use.		
Sterilization method	Ethylene oxide (ETO).		
Supply	Sterile, single container. Device which includes a bulb syringe, latex free, that facilitates the positioning and extraction of the ocular protector.		
Shelf-life	4.5 years		

* Product with CE marking



REV: 05. 07/2022