



OCULFIT LÁMINA

Suelo Orbital

Nombre de marca	OCULFIT		
Referencias	8438-POI	7210-POI	6330-POI
Medidas (LxAxH)	50 x 38 x 0,4 mm ULTRAFINO	50 x 38 x 0,85 mm FINO	50 x 38 x 1,5 mm ESTÁNDAR
Descripción e indicaciones	La lámina de suelo orbital OCULFIT está fabricada en polietileno poroso de alta densidad. Su diseño de estructura de poros abiertos interconectados permite la vascularización y el crecimiento de tejidos dentro de él. Está demostrado que se trata de un material estable y biocompatible con los tejidos, lo que proporciona a los cirujanos un implante seguro para la reconstrucción orbitaria. Este modelo de implante está indicado en: - Reconstrucción de la órbita para prevenir o tratar la diplopía causada por la disminución del volumen orbital o por un atrapamiento del músculo recto inferior. - Reconstrucción del suelo orbital evitando los problemas funcionales del enoftalmos, que provoca un aumento en volumen orbital.		
Material	Polietileno poroso de alta densidad (HDPE).		
Método de esterilización	Óxido de etileno (ETO).		
Presentación	Envase unitario en doble bolsa Tyvek, estéril.		
Caducidad	5 años		

* Productos con marcado CE



REV: 05. 07/2022



OCULFIT SHEET

Orbital Implant

Brand name	OCULFIT		
References	8438-POI	7210-POI	6330-POI
Measures (LxAxH)	50 x 38 x 0,4 mm ULTRATHIN	50 x 38 x 0,85 mm THIN	50 x 38 x 1,5 mm STANDARD
Description and Indications	The OCULFIT orbital implant is made of high-density polyethylene. Its interconnected open pore structure design allows vascularisation and tissue growth inside it. It is proven to be a stable material that is biocompatible with tissues, providing surgeons with a safe implant for orbital reconstruction. This implant model is suitable for: - Reconstructing the orbit to prevent or treat diplopia caused by decreased orbital volume or by an entrapment of the inferior rectus muscle. - Reconstructing the orbital floor avoiding the functional problems of the enophthalmos, which causes an increase in orbital volume.		
Material	High Density Poly Ethylene (HDPE).		
Sterilization method	Ethylene Oxide (ETO).		
Supply	Unit packaged in dual barrier Tyvek bag, sterile.		
Shelf-life	5 years		

* Product with CE marking



REV: 05. 07/2022