




ANILLO DE SIMBLEFARON


Oculoplastia

Nombre de marca	AJL ANILLO DE SIMBLÉFARON		
Referencias	SF-S-AJL	SF-M-AJL	SF-L-AJL
Tamaño	Pequeño: 22 x 19 x 1.2 mm	Mediano: 24 x 21 x 1.2 mm	Grande: 25 x 22 x 1.2 mm
Descripción e Indicaciones	El ANILLO DE SIMBLÉFARON AJL es un dispositivo médico diseñado con una geometría específicamente adaptada al globo ocular. Se utiliza para evitar la adhesión total o parcial entre la cara interna del párpado y el globo ocular, así como para el mantenimiento del fondo de saco conjuntival.		
Material	Acrílico de alta resistencia y acabado superficial apto para el uso ocular.		
Método de esterilización	Óxido de etileno (ETO).		
Presentación	Envase unitario, estéril.		
Caducidad	4,5 años		

* Producto con marcado CE



REV: 05. 07/2022



SYMBLEPHARON RING

Oculoplasty

Brand name	AJL SYMBLEPHARON RING		
References	SF-S-AJL	SF-M-AJL	SF-L-AJL
Size	Small: 22 x 19 x 1.2 mm	Medium: 24 x 21 x 1.2 mm	Large: 25 x 22 x 1.2 mm
Description and Indications	AJL SYMBLEPHARON RING is a medical device designed with a geometry specifically adapted to the eyeball. It's used to avoid total or partial adhesion between the internal face of the eyelid and the eyeball, as well as maintaining the depth of the conjunctival bag.		
Material	High resistance acrylic material suitable for ocular use.		
Sterilization method	Ethylene Oxide (ETO).		
Supply	Sterile, single container.		
Shelf-life	4.5 years		

* Product with CE marking



REV: 05. 07/2022