

AS – IOL


Lente Intraocular

Referencia	LLASY60
Nombre de marca	AS-IOL
Descripción	Lente intraocular acrílica asférica. Diseño monobloque. Cámara posterior.
Material	Acrílico, hidrofílico.
Filtro	Filtro UV / Filtro natural para la luz azul.
Diámetro de la óptica	6.0 mm
Longitud total	10.75 mm
Angulación	0°
Rango de dioptrías	-5 D a +40 D en incrementos de 1 D / +12 D a +25 D en incrementos de 0,5 D
Borde	Cuadrado (360°)
Contenido de agua	25.5%
Constante A / ACD	Constante A (óptica): 118 Constante A (U/S): 117,3 ACD: 4,56
Método de esterilización	Vapor.
Presentación	Envase unitario, estéril
Caducidad	3 años

* Producto con marcado CE



REV: 05. 07/2022



AS – IOL

Lente Intraocular

Reference	LLASY60
Brand name	AS-IOL
Description	Aspherical acrylic intraocular lens. One piece design. Posterior chamber.
Material	Hydrophilic acrylic.
Filter	UV Filter / natural blue light filter.
Optical diameter	6.0 mm
Total length	10.75 mm
Angle	0°
Diopters range	-5 D to +40 D in steps of 1 D / +12 D to +25 D in steps of 0.5 D
Edge	Square (360°)
Water content	25.5%
Constant A / ACD	Constant A (optical): 118 Constant A (U/S): 117.3 ACD: 4.56
Sterilization method	Steam.
Supply	Sterile, single unit
Shelf-life	3 years

* Product with CE marking



REV: 05. 07/2022