



AJL Ophthalmic, S.A. is the only Spanish company (July 1992) dedicated to the R&D, manufacture and marketing of medical devices in the field of Ophthalmology.

The company has reached a leading market position thanks to the fulfillment of very exigent Regulatory and Quality Control Requirements, 3000 m2 facilities with cutting-edge technology as well as the most qualified personnel for the manufacture of its products.

AJL has an extensive sales network with presence in more than 80 countries.



AJL VET



CANINE IOLS

10.25MM TXY651025
11.25MM TXY651125
12.25MM TXY651225
13.25MM TXY651325
14.25MM TXY651425
15.25MM TXY651525

41D 6.5 OPT
Disposable injector is supplied



FELINE IOLS

13MM CAT701300
53.5D 7.0 OPT

Disposable injector is supplied



PMMA CANINE IOLS

14MM - 6.5 OPT - P651400
16MM - 7.0 OPT - P701600

41D



TRYPAN BLUE

0.6mg/ml BLUEVET

5x 1ML SYRINGES/BOX

Local regulations may apply



CTRS

11MM ACV001102
13MM ACV001302
15MM ACV001502
13MM* ACV001302PL
16MM* ACV001602PL

*PRE-LOADED



VISCOELASTICS

NaHA 1.4% VV14
NaHA 2% VV20
NaHA 3% VV30
HPMC 2% VCELL

5x 2ML SYRINGES/BOX



SILICONE OIL

5000 cSt S5000VET

10ML SYRINGE



MEMBRANES

2x2cm 100um MCV22100

2x2cm 200um MCV22200

COLLAGEN TYPE I



AJL Ophthalmic S.A.
Parque Tecnológico de Álava
C/ Ferdinand Zeppelin 1
01510 Miñano, Álava - Spain

T. +34 945 298 256

F. +34 945 298 209

ajlsa@ajlsa.com

www.ajlsa.com



AJL Ophthalmic S.A.
Parque Tecnológico de Álava
C/ Ferdinand Zeppelin 1
01510 Miñano, Álava - Spain

T. +34 945 298 256

F. +34 945 298 209

ajlsa@ajlsa.com

www.ajlsa.com



The future of vision

AJL PRODUCT PORTFOLIO

CORNEA




Keratoconus Intrastromal Corneal Rings



Intacs[®] SK (6mm optical zone)



AFR (5 mm optical zone)
AFR6 (6 mm optical zone)



AJL pro+

CW (clockwise/5 mm optical zone)
CW6 (clockwise/6 mm optical zone)
CCW (counterclockwise/5 mm optical zone)
CCW6 (counterclockwise/6 mm optical zone)

VACUUM SYSTEM

KV2000



Medical device designed to stabilize the cornea and facilitate the tunneling process during the intracorneal rings implantation surgery



BRIEFCASE INCLUDED

Dimensions: 46x33x15
Weight: 5,10 Kg.

Endo-keratoprothesis



ENDO-K PRO

EQP75 7.50mm
EQP80 8.00mm
EQP85 8.50mm
EQP90 9.00mm

Pre-descemetic Corneal Ring



NEORING

AQ080 8.00mm
AQ085 8.50mm
AQ090 9.00mm

CATARACT



Iris expander



CANABRAVA RING
CANA-RING

PMMA Intraocular lenses



AJL PMMA
P651300
(Posterior chamber)



AJL PMMA
A601250
(Anterior chamber)

Capsule stain solution



AJL BLUE
0,06%

Hydrophilic Intraocular lenses



AIALA
F601250



AS-IOL
LLASY60

Hydrophobic Intraocular lenses



AIALA DRY
LLASHP60



AIALA DRY PRELODED
LLASHP60-PL

CAPSULAR TENSION RING



AC
ACS
AC-PL
ACS-PL

Open diameter	10mm	11mm	12mm	12mm	13mm	14mm
Close diameter	9mm	10mm	10mm	10,5mm	11mm	12mm

Viscoelastic solutions



AJL VISC

- AJL VISC 1%
- AJL VISC 1,4%
- AJL VISC 2%
- AJL VISC 3%
- AJL CELL 2%

GLAUCOMA



Scleral implant



ESNOPER
V2000

Uveoscleral implant



ESNOPER CLIP
VCLIP

NPDS surgical instruments (non-penetrating deep sclerectomy)



NPDS KIT

RETINA



Macular Buckle



NPB

ILM stain solution



AJL BBG

Silicone oil



SIOBAL
S1000
S2000
S5000

SIOBAL

OCULOPLASTY



Orbital Implants



OCULFIT SPHERICAL

6316-POI 14mm	6317-POI 20mm
6326-POI 16mm	6703-POI 21mm
6327-POI 18mm	6322-POI 22mm
8702-POI 19mm	8704-POI 23mm



OCULFIT CONICAL

80016-POI 16mm
80018-POI 18mm
80020-POI 20mm



OCULFIT SHEET

8438-POI ULTRA THIN
7210-POI THIN
6330-POI STANDARD



OCULAR SHIELD

PROCUL-S-AJL
PROCUL-M-AJL
PROCUL-L-AJL



SYMBLEPHARON RING

SF-S-AJL
SF-M-AJL
SF-L-AJL



PERFORATED COMFORMER

CF-S-AJL
CF-M-AJL
CF-L-AJL



ALOS PALPEBRAL IMPLANTS

GOLD TH (0,6 a 1,8)
PLATINUM PHT (0,6 a 1,8)