



ALOS



Palpebral Implants

Brand name	ALOS
Models	GOLD (TH - G) and PLATINUM (PTH - P)
Description and Indications	• Biocompatible materials • Due to their fine, spherical design, the implants adapt to the curvature of the eyeball to achieve a better aesthetic result. • Rounded edges to avoid the extrusion of the implant in the long-term. • Suture holes with fixation channels. • Purity of 99,9%(Platinum) and 99,9%(Gold). • Sterile. <u>Indication:</u> Eyelid implants are designed for the correction of lagophthalmos.



Brand name	ALOS SIZER SET
References	SET MEDIDOR PESAS TH SET MEDIDOR PESAS G
Weight (g)	0.6-0.8-1.0-1.2-1.4-1.6

GOLD

References	TH 06	TH 08	TH 10	TH 12	TH 14	TH 16	TH 18	G 06	G 08	G 10	G 12	G 14	G 16	G 18	G 20	G 22	G 24	G 26
Weight (g)	0.6	0.8	1.0	1.2	1.4	1.6	1.8	0.6	0.8	1.0	1.2	1.4	1.6	1.8	2.0	2.2	2.4	2.6
Length (mm)	10.5	12.5	15.5	17.5	19.5	21.5	24.5	7.8	10.3	12.3	13.9	15.4	16.9	18.9	20.5	20.8	22.3	23
Width (mm)	5.7	5.7	5.7	5.7	5.7	5.7	5.7	5	5	5	5	5	5	5	5	5	5	5
Thickness (mm)	0.6	0.6	0.6	0.6	0.6	0.6	0.6	1.1	1.1	1.1	1.1	1.1	1.1	1.1	1.1	1.1	1.1	1.1
Sterilization method	Steam																	
Shelf-life	3 years																	

PLATINUM

References	PTH06	PTH08	PTH 10	PTH 12	PTH 14	PTH 16	PTH 18	P 06	P 08	P 10	P 12	P 14	P 16	P 18
Weight (g)	0.6	0.8	1.0	1.2	1.4	1.6	1.8	0.6	0.8	1.0	1.2	1.4	1.6	1.8
Length (mm)	10.2	12.4	15.0	17.0	18.8	21.0	24.0	7.6	10	11.9	13.5	15.0	16.3	17.4
Width (mm)	5.7	5.7	5.7	5.7	5.7	5.7	5.7	5	5	5	5	5	5	5
Thickness (mm)	0.6	0.6	0.6	0.6	0.6	0.6	0.6	1.1	1.1	1.1	1.1	1.1	1.1	1.1
Sterilization method	Steam													
Shelf-life	3 years													

* Product with CE marking



REV: 05. 07/2022