



AS – IOL

Lente Intraocular

Reference	LLASY60
Brand name	AS-IOL
Description	Aspherical acrylic intraocular lens. One piece design. Posterior chamber.
Material	Hydrophilic acrylic.
Filter	UV Filter / natural blue light filter.
Optical diameter	6.0 mm
Total length	10.75 mm
Angle	0°
Diopters range	-5 D to +40 D in steps of 1 D / +12 D to +25 D in steps of 0.5 D
Edge	Square (360°)
Water content	25.5%
Constant A / ACD	Constant A (optical): 118 Constant A (U/S): 117.3 ACD: 4.56
Sterilization method	Steam.
Supply	Sterile, single unit
Shelf-life	3 years

* Product with CE marking



REV: 05. 07/2022