



AIALA DRY

Intraocular Lenses

Brand name	AIALA DRY	AIALA DRY PRELOADED
References	LLASHP60	LLASHP60-PL
Description	Monophocal, aspheric, posterior chamber intraocular lens.	
Design	Single piece. C-Loop with antistick design. Square Edge (360°)	
Material	Hydrophobic acrylic, UV and blue light filter. Glistening-Free.	
Optical diameter	6.00 mm	
Total length	13.00 mm	
Angle	0°	
Diopters range	-5 D to +40 D (in steps of 1.0 D) / +10 D to +30 D (in steps of 0.5 D)	
Water content	<4%	
Constant A / ACD	A-constant (optical): 119.7 A-constant (U/S): 119.2 ACD: 5.66	
Sterilization method	Ethylene Oxide (ETO).	
Supply	Single unit, sterile.	
Shelf-life	3 years	

* Product with CE marking



* Preloaded intraocular lens injector
Ref.: LLASHP60-PL



REV: 05. 07/2022