



AIALA



IOLs

Reference	F601250
Brand name	AIALA
Description	Hydrophilic intraocular lens. One piece design. Posterior chamber.
Desing	C Haptics.
Material	Hydrophilic acrylic.
Filter	UV-A
Optical diameter	6.0 mm
Total length	12.50 mm
Angle	0°
Diopters range	-6 D to +40 D in steps of 1 D / +10 D to +25 D in steps of 0.5 D
Edge	Square (360°)
Water content	25.5%
Constant A / ACD	Constant A: (optical) 118 Constant A: (u/s) 117.3 ACD: 4.9
Sterilization method	Steam.
Supply	Single unit, sterile.
Injector system recommended	Medicel.
Shelf-life	3 years.

* Product with CE marking



REV: 05. 07/2022