



PERFORATED CONFORMER

Perforated Conformer

<i>Brand name</i>	AJL PERFORATED CONFORMER		
<i>References</i>	CF-S-AJL	CF-M-AJL	CF-L-AJL
<i>Size</i>	Small: 22 x 19 x 1.2 mm	Medium: 24 x 21 x 1.2 mm	Large: 25 x 22 x 1.2 mm
<i>Description and Indications</i>	<p>The AJL PERFORATED CONFORMER is a medical device designed to avoid adherence between the tarsal and bulbar conjunctiva. It is transparent and has two orifices for ventilation at the centre.</p> <p>AJL PERFORATED CONFORMER is used to achieve the separation between tarsal and bulbar conjunctiva. The device prevents the shrinkage of the cavity and maintains the orbital space properly after enucleation/evisceration.</p>		
<i>Material</i>	High resistance acrylic material suitable for ocular use.		
<i>Sterilization method</i>	Ethylene Oxide (ETO).		
<i>Supply</i>	Sterile, single container.		
<i>Shelf-life</i>	4.5 years		

* Product with CE marking



REV: 05. 07/2022