



INTACS® SK

Intrastromal Corneal Implants

<i>Model</i>	Intacs®SK (optical area of 6 mm).
<i>Brand name</i>	Intacs®SK
<i>Definition</i>	Intrastromal corneal implants
<i>Indications</i>	Consult nomogram for each case. Customized advice: keratoconus@ajlsa.com
<i>Product description</i>	Intacs®SK are intrastromal corneal segments designed for its implantation the corneal stroma through a radial small incision. The edges are blunt in order to reduce halos and other undesirable visual effects.
<i>Material</i>	Poly-methylmethacrylate (PMMA)
<i>Sterilization method</i>	Ethylene oxyde (ETO)
<i>Supply</i>	Single sterile packaging
<i>Shelf-life</i>	5 years

INTACS® SK (6.0 mm) CORNEAL IMPLANT

THICKNESS-ARC	REFERENCE NUMBER
0,250mm-90º	12300-250-90
0,300mm-90º	12300-300-90
0,350mm-90º	12300-350-90
0,400mm-90º	12300-400-90
0,450mm-90º	12300-450-90

0,300mm-130º	12300-300-130
0,350mm-130º	12300-350-130
0,400mm-130º	12300-400-130
0,450mm-130º	12300-450-130

0,210mm-150º	12300-210-150
0,250mm-150º	12300-250-150
0,300mm-150º	12300-300-150
0,350mm-150º	12300-350-150
0,400mm-150º	12300-400-150
0,450mm-150º	12300-450-150
0,500mm-150º	12300-500-150

* OTHER THICKNESS AVAILABLE UNDER PETITION.



REV: 05. 07/2022