




AJL BBG



Solución Oftalmológica

Nombre de marca	AJL BBG
Referencia	AJL BBG 0,025%
Descripción e indicaciones	AJL BBG 0,025% es una solución oftalmológica estéril de azul brillante G, no pirógena, isotónica y preparadas para su uso en la cirugía vitreoretinal. La solución BBG es capaz de teñir eficazmente la membrana limitante interna (MLI). <u>Indicaciones:</u> AJL BBG está indicado en patologías maculares como en ojos con agujero macular, edema macular o retinopatía diabética. <u>Especificaciones:</u> <ul style="list-style-type: none">• AJL BBG está destinado a ser aplicado en la cámara vítreo diferenciando así la membrana limitante interna (MLI) del tejido subyacente de la retina y de la membrana epirretiniana. El colorante no penetra en la retina, permitiendo la visualización de la MLI.• El exceso de colorante se elimina mediante técnicas de aspiración/irrigación con solución salina.
Composición	Solución oftalmológica estéril de azul brillante G al 0,025%.
Método de esterilización	Vapor.
Presentación	Jeringa de vidrio precargada de 1ml, estéril, de un solo uso con conexión Luer Lock.
Caducidad	3 años


€ 2797



V-04-01-AAT-IMP REV.00



AJL BBG



Ophthalmic Solution

<i>Brand name</i>	AJL BBG
<i>Reference</i>	AJL BBG 0.025%
<i>Description and Indications</i>	<p>AJL BBG 0.025% is an ophthalmic sterile, non-pyrogenic, isotonic solution, ready for use in vitreoretinal surgery. The BBG solution is capable of effectively staining the inner limiting membrane (ILM).</p> <p><u>Indications:</u></p> <p>AJL BBG is indicated for macular pathologies such as eyes with macular hole, macular edema or diabetic retinopathy.</p> <p><u>Specifications:</u></p> <ul style="list-style-type: none">• AJL BBG is intended to be applied in the vitreous chamber thereby differentiating the internal limiting membrane (ILM) from the underlying retinal tissue and the epiretinal membrane. The dye does not penetrate the retina, allowing visualization of the ILM.• Excess dye is removed by aspiration irrigation techniques with saline solution.
<i>Composition</i>	Brilliant blue G sterile ophthalmic solution at 0.025%
<i>Sterilization method</i>	Steam.
<i>Supply</i>	Single-use, sterile, 1ml pre-filled glass syringe with Luer Lock connection.
<i>Shelf-life</i>	3 years

€ 2797



V-04-01-AAT-IMP REV.00