




AJL BBG



Ophthalmic Solution

Brand name	AJL BBG
Reference	AJL BBG 0.025%
Description and Indications	<p>AJL BBG 0.025% is an ophthalmic sterile, non-pyrogenic, isotonic solution, ready for use in vitreoretinal surgery. The BBG solution is capable of effectively staining the inner limiting membrane (ILM).</p> <p><u>Indications:</u></p> <p>AJL BBG is indicated for macular pathologies such as eyes with macular hole, macular edema or diabetic retinopathy.</p> <p><u>Specifications:</u></p> <ul style="list-style-type: none">• AJL BBG is intended to be applied in the vitreous chamber thereby differentiating the internal limiting membrane (ILM) from the underlying retinal tissue and the epiretinal membrane. The dye does not penetrate the retina, allowing visualization of the ILM.• Excess dye is removed by aspiration irrigation techniques with saline solution.
Composition	Brilliant blue G sterile ophthalmic solution at 0.025%
Sterilization method	Steam.
Supply	Single-use, sterile, 1ml pre-filled glass syringe with Luer Lock connection.
Shelf-life	3 years

€ 2797

