



AJL BLUE

Ophthalmic Solution

Brand name	AJL BLUE
Reference	AJL BLUE 0.06%
Description and indications	<p>AJL BLUE is an ophthalmic solution of Trypan Blue at 0.06%, capable of effectively staining the anterior lens capsule, facilitating the capsulorhexis process in cataract surgery. It is a buffered, isotonic, sterile and non-pyrogenic ophthalmic stain with physiological pH.</p> <p>AJL BLUE is indicated for eyes with dense, mature cataracts, poor fundus reflex or small pupils.</p> <ul style="list-style-type: none">• AJL BLUE is intended to be applied in the anterior chamber of the eye, on the anterior capsule of the lens, thus staining any portion of the capsule that comes into contact with the dye.• Excess dye is removed by aspiration irrigation techniques with saline solution.
Composition	Trypan blue sterile ophthalmic solution at 0.06%
Sterilization method	Steam
Supply	Single-use, sterile, 1ml pre-filled glass syringe with Luer Lock connection and 27G sterile cannula.
Purity content	>98%
Shelf-life	3 years

CE 2797

