




AJL PMMA

PMMA Intraocular Lenses

Reference	P651300 	A601250 
Brand name	AJL-PMMA	AJL-PMMA
Model	Posterior chamber, spherical, one piece.	Anterior chamber, spherical, one piece.
Material	Polymethyl Methacrylate (PMMA)	Polymethyl Methacrylate (PMMA)
Filter	UV	UV
Optical diameter	6,5 mm	6,0 mm
Total length	13 mm	12.5 mm
Angle	10°	16°
Diopters range	-10,0 D to +40,0 D in steps of 1,0 D / +12,0 D to +25,0 D in steps of 0,5 D	+6,0 D to +30,0 D in steps of 1,0 D
Edge	Square (360°)	Square (360°)
Supply	Sterile, single unit.	Sterile, single unit.
Constant A / ACD	Constant A (optical): 118.6 Constant A (U/S): 118 ACD: 4.96	Constant A (optical): 116.2 Constant A (U/S): 115.5 ACD: 3.51
Sterilization method	Ethylene Oxide (ETO).	Ethylene Oxide (ETO).
Shelf-life	4.5 years	4.5 years

CE 0318



V-04-01-AAJ-IMP REV.00