



AS – IOL

Lente Intraocular

Reference	LLASY60
Brand name	AS-IOL
Description	Acrylic posterior chamber monofocal monobloc monofocal intraocular lens
Haptic desing	Closed loop haptics with fenestrated plate shape
Material	Hydrophilic acrylic.
Filter	UV Filter / natural blue light filter.
Optical diameter	6,00 mm
Total length	10,75 mm
Angle	0°
Diopters range	-5,0 D to +40,0 D in steps of 1,0 D / +12,0 D to +25,0 D in steps of 0,5 D
Edge	Square (360°)
Water content	25.5%
Constant A / ACD	Constant A: (optical) 118 Constant A: (u/s) 117.30 ACD: 4.56
Sterilization method	Steam
Supply	Single unit, sterile.
Injector system recommended	Viscoject 2.2
Shelf-life	3 years.

€ 2797



V-04-01-AAM-IMP REV.00