



AIALA DRY

Intraocular Lenses

| | | |
|-----------------------------|--|----------------------------|
| <i>Brand name</i> | AIALA DRY | AIALA DRY PRELOADED |
| <i>References</i> | LLASHP60 | LLASHP60-PL |
| <i>Description</i> | Single Piece, monophocal, posterior chamber intraocular lens | |
| <i>Haptic desing</i> | C-Loop with antistick design | |
| <i>Material</i> | Hydrophobic acrylic, UV and blue light filter. Glistening-Free. | |
| <i>Optical diameter</i> | 6,00 mm | |
| <i>Total length</i> | 13,00 mm | |
| <i>Angle</i> | 0° | |
| <i>Diopters range</i> | -5,00 D to +40,00 D (in steps of 1,00 D) / +10,00 D to +30,00 D (in steps of 0,50 D) | |
| <i>Water content</i> | <4% | |
| <i>Constant A / ACD</i> | A-constant (optical): 119,7 A-constant (U/S): 119,2 ACD: 5,66 | |
| <i>Sterilization method</i> | Ethylene Oxide (ETO). | |
| <i>recommended injector</i> | Accuject 2.2 (For LLASHP60) | |
| <i>Supply</i> | Single unit, sterile. | |
| <i>Shelf-life</i> | 3 years | |



* Preloaded intraocular lens injector
Ref.: LLASHP60-PL

CE 2797



V-04-01-AAN-IMP REV.00