



PERFORATED CONFORMER

Perforated Conformer

Brand name	AJL PERFORATED CONFORMER		
References	CF-S-AJL	CF-M-AJL	CF-L-AJL
Size	Small: 22 x 19 x 1.2 mm	Medium: 24 x 21 x 1.2 mm	Large: 25 x 22 x 1.2 mm
Description and Indications	<p>The Perforated Conformer is an oculoplastic product designed to achieve the separation between the tarsal conjunctiva (covers the inner side of the eyelid) and bulbar conjunctiva (covers the eyeball), which is the moist tissue that covers the eye in its outermost part, and the eyelid in its inner part, with the intention of avoiding adherence between one and the other, and at the same time preserving the space preventing cavitary contraction.</p> <p>Its use is indicated as a complementary product after evisceration or enucleation surgery, from the moment the surgery is completed until the patient has his definitive prosthesis (as long as this period does not exceed 30 days); thus avoiding closure or adhesions during the healing process.</p>		
Material	It is made of a highly resistant thermoplastic translucent material of polyethylene and copolyester (PTEG) and biocompatible.		
Sterilization method	Ethylene Oxide (ETO).		
Supply	Sterile, single container.		
Shelf-life	4.5 years		

CE 2797

