



SYMBLEPHARON RING



Oculoplasty

Brand name	AJL SYMBLEPHARON RING		
References	SF-S-AJL	SF-M-AJL	SF-L-AJL
Size	Small: 22 x 19	Medium: 24 x 21	Large: 25 x 22
Description and Indications	<p>The Simblefaron AJL Ring is an oculoplastic product designed to prevent conjunctival fornix retraction. The central width of the hole prevents corneal contact while preserving the cornea.</p> <p>The Simblefaron ring is designed to prevent conjunctival fornix retraction while preserving the cornea. It is placed under the patient's eyelid in order to prevent the appearance of what is called "symblepharon" (total or partial adhesion between the inner face of the eyelid and the eyeball as a consequence of any scarring process). This anomaly can affect both eyelids, although it is more common in the lower eyelid.</p>		
Material	It is made of a highly resistant thermoplastic translucent material of polyethylene and copolyester (PTEG) and biocompatible.		
Sterilization method	Ethylene Oxide (ETO).		
Supply	Sterile, single container.		
Shelf-life	4.5 years		

€ 2797



V-04-01-AAZ-IMP REV.00