



SIOBAL






Silicone Oil

Brand name	SIOBAL		
Reference	S1000	S2000	S5000
Description and indications	SIOBAL silicone oil is a viscous, transparent, homogeneous, sterile and non-pyrogenic liquid, which is made up of high purity polydimethylsiloxane. Tamponade agent for retinal detachment. Is used as a tamponade product in vitreoretinal surgeries to treat long-term complex retinal detachments, especially in cases of: <ul style="list-style-type: none">• Retinal detachment• Proliferative vitreoretinopathy.• Perforating wounds.• Proliferative diabetic retinopathy.		
Viscosity (1s ⁻¹ , 25°C)	950 - 1.500 mPa·s	1.500 - 2.500mPa·s	4.500 - 5.500 mPa·s
Composition	100% silicone oil (polydimethylsiloxane)		
Refractive index	(25°C) 1,40		
Purity	>95%		
Sterilization method	Steam.		
Supply	Single-use, sterile, 10ml pre-filled glass syringe with Luer Lock connection		
Shelf-life	3 years		

CE 0318








SIOBAL PACK

Brand name	SIOBAL PACK: SIOBAL Silicone Oil + Inyection Tubing Set (Sterile)											
References	PS1000 *				PS2000 *				PS5000 *			
Models	A	OA	OR	B	A	OA	OR	B	A	OA	OR	B
	*Selection of compatible connection system for the following equipments: (A) Alcon: Accurus, Constellation (B) Bausch&Lomb: Millenium, Stellaris											
Description and indications	SIOBAL silicone oil injection /extraction system. This system in connected to vitrectomy equipment so as to facilitate the SIOBAL silicone oil injection/extraction process.											
Viscosity- oil	950 - 1.500 mPa.s				1.500 - 2.500mPa.s				4.500 - 5.500 mPa.s			
Content	Cannula 20G		Cannula 23G		Cannula 25G		Connection system for injection / extraction		10ml oil prefilled syringe			
												
Supply	5 packs per box. Products sterilized and packed in their original containers by the manufacturer.											
Quality assurance	CE 0318											



SIOBAL SET

Brand name	SIOBAL SET: Silicone Oil Injection / Extraction Tubing Set (Sterile)				
Reference	SIE*				
Models	A	OA	OR	B	
	*Selection of compatible connection system for the following equipments: (A) Alcon: Accurus, Constellation Optikon: (OA) Assistant, (OR) R-Evolution. (B) Bausch&Lomb: Millenium, Stellaris				
Content	Cannula 20G	Cannula 23G	Cannula 25G	Connection system for injection / extraction	Syringe 10ml
				 SIE-A SIE-OA SIE-OR SIE-B	
Supply	5 sets per box. Products sterilized and packed in their original containers by the manufacturer.				
Quality assurance	CE 0318				