



NPB



Macular Buckle

Brand name	NPB
Definition	Macular Buckle
Indications	<ul style="list-style-type: none">• MYOPIC TRACTION MACULOPATHY in stage of maculoschisis and macular detachment, with and without macular hole
Product description	The NPB macular buckle is made of medical grade PMMA and coated with a thin biocompatible silicon layer. It has an indented surface with a spherical shell at the top and an arm which length varies depending on the axial lenght of the patient's eye. The silicone belt inserted in the arm of the implant facilitates the fixation of a light probe, through its reflection, allows the correct positioning in the desired indentation site
Total length of the arm	The dimensions of the implant will depend on the axial length of the patient's eye*
Material	Polymethyl Methacrylate (PMMA) coated with silicone of medical grade
Sterilization method	Ethylene Oxide (ETO)
Supply	Sterile, single unit.
Shelf-life	4.5 years

* To order the specific NPB, please contact AJL Ophthalmic or the distributor in your region.

