



AJL OOKP

Osteo-odonto-keratoprosthesis

<i>Brand name</i>	AJL OOKP
<i>Reference</i>	AJL-OOKP
<i>Product description and indications</i>	<p>The osteo-odonto-keratoprosthesis is a medical device implantable in the cornea with a transparent cylinder shape supported by an autologus tissue (canine tooth or tibia bone). It enables to restore sight to an eye in cases of corneal opacity that avoids vision and a new corneal transplant is not possible:</p> <ul style="list-style-type: none"> - Autoimmune diseases such as Stevens-Johnson syndrome and ocular pemphigoid - Severe dry eye - Severe chemical burns affecting both eyes
<i>Implant specifications</i>	Custom made device
<i>Ordering documentation</i>	<ul style="list-style-type: none"> - Complete the order form specifying your patient's axial length and additional device data. - Anterior segment OCT
<i>Material</i>	Polymethyl Methacrylate (PMMA)
<i>Sterilization method</i>	Ethylene Oxide (ETO)
<i>Supply</i>	Single use packaging
<i>Shelf-life</i>	4.5 years

*To order the specific Osteo-odonto-Keratoprosthesis, please contact AJL Ophthalmic or the distributor in your region.

