




VAMA



Intubation Cannula

Brand name	VAMA
Reference	VMA
Product description and indications	The VAMA Cannula is indicated as a guide for fibroscopic intubation via the oral route with endotracheal tubes with a maximum internal diameter of 8.0 mm. It can be used to intubate awake or anaesthetised patients, with or without muscle relaxants.
Product specifications	<ul style="list-style-type: none">-Anatomical shape.-It has a laser mark on the midline of the cannula body to facilitate insertion of the fibroscope, allowing easy guidance of itself.-Its sliding lid facilitates safe removal of the cannula once the tube is in its functional intratracheal position.
Dimensions (length x inside diameter x height)	105.74 x 16.50 x 42.00 mm
Material	ABS
Supply	Single-use medical device, latex-free and non-sterile. Unitary packaging

CE



V-04-01-ABU-IMP REV.01