

TF series

ophthalmic trial frames

TF-1 (prev. TF-2 Pediatric)

Trial frame with special, adjustable nose pads perfect for examination of children. It's characterised by high durability & excellent quality.




TF-2

Light-weight trial frame designed to be used with both, adults and children (two different size nose rests included), assures full adjustment against nose and ears position.




TF-4

Solid trial frame designed to be used with adults, assures full adjustment against nose and ears position.



TF-6

An essential tool for optometrists and ophthalmologists, designed to facilitate precise refractive assessments. Crafted with durable yet lightweight materials, ensures comfort for patients during extensive fittings.



Technical specification	TF-1	TF-2	TF-4	TF-6
Weight ($\pm 5g$)	47g	61g	80g	55g
PD adjustment range (pupil spacing)	48 - 60mm	50 - 80mm	48 - 80mm	50 - 80mm
Step of the axis scale	5°	5°	5°	5°
Range of lens rotation	0 - 360°	0 - 360°	0 - 360°	0 - 360°
Temple tips angle adjustment range	$\pm 15^\circ$	$\pm 10^\circ$	$\pm 10^\circ$	$\pm 10^\circ$
Vertex distance adjustment range	n/a	0 - 13mm	0 - 10mm	0 - 13mm
Nose pad adjustment range	0 - 4mm	0 - 23mm	0 - 14mm	0 - 23mm
Max. number of trial lenses	4	5	4	5
Packaging	1box; 19x11x6cm, 0.25kg			