

SC-1

surgical table



TECH SPEC

- Full range of patient's positioning (from sitting to lying, incl. Trendelenburg position) what gives unprecedented opportunity to position the patients at any age in a comfortable way
- A unique design of the headrest gives the surgeon precise control over the patient's head during surgical procedures
- Detachable armrests and cushioning supports, as well as stainless steel push handles
- Ease of movement around the operation theatre
- Central locking gear which assures stability during the surgery AC/DC power supply operated
- Battery pack - AC/DC power supply operated
- 4 universal guide rails (EU)
- A selection of non-stitched upholstery available without a surcharge

Technical specification

Height adjustment range	600 - 750mm
Trendelenburg position	ca. 15°
Max. load capacity	150kg
Supply voltage*	~100 - 240V AC, 50/60Hz
Battery pack	24V, 12Ah
Battery charging	24V ±0.5V, 1.2A
Working time	up to 40h
Packaging	1 pallet; 120x80x160cm; 164kg (gross)

Dimensions

Length (flattened)	1 740 - 1 840mm
Length (sitting position)	1 120mm
Width	874mm
Height	1 240 - 1 390mm

*- The power supply can support both 100VAC and 240VAC power sources. Voltage value detection build-in.



vanilla



mint



juniper



very berry