

Electric Lift Tables

TT2C-800/TT2C-1000

The TT2C-800 and TT2C-1000 are vertically adjustable electric tables with two columns. Thanks to their size and stable stand they are ideal for placing two eye diagnostic devices. Also, the access for wheelchair users is quite convenient. The maximum load of the tables is 90 kg. Optional shelves for cables and accessories are available.

TTUD-1000

This model is designed for our ultrasound line. It is prepared with a drawer for storing the probes and a holder for wipes and cleaning liquids. The maximum load is 65 kg.

Specifications

	TT2C-800	TT2C-1000	TTUD-1000
DIMENSIONS			
Max. tabletop load	90 kg	90 kg	65 kg
Vertical movement range	300 mm	300 mm	302 mm
Tabletop dimensions	1000 x 500 mm	1240 x 550 mm	580 x 545 mm
Min. / max. tabletop height	608-908 mm	608-908 mm	702-1004 mm
Weight	33.6 kg	42 kg	34.6 kg

TT2C-800



TT2C-1000

Excellent features

TT2C-800 / TT2C-1000

- + Tabletops:
TT2C-800: 1000 (W) x 500 (D) mm
TT2C-1000: 1240 (W) x 550 (D) mm
- + Vertical adjustment: 608 to 908 mm
- + Maximum load: 90 kg
- + Two-column design

TTUD-1000

- + Tabletop: 580 (W) x 545 (D) mm
- + Vertical adjustment: 702 to 1004 mm
- + Maximum load: 65 kg
- + Fits perfectly for ultrasound line

TTUD-1000



TOMEY EUROPE
TOMEY GMBH

tomey.de

TOMEY GmbH is the European headquarter of TOMEY Corporation,
2-11-33 Noritakeshinmachi Nishi-Ku, Nagoya, 451-0051, Japan

2023/08 – subject to change without notice

Always read and follow the instructions for use.
Not all products, services or offers are approved or offered in every market. Please note that the current status of approval for the labelling, instructions and contents of the brochure may vary from one country to another.

