

OC-1

ophthalmic operator's chair



TECH SPEC

- Maximum comfort for the ophthalmic surgeons during cataract surgery and other procedures ensured
- Electrical height adjustment thanks to ergonomically positioned foot switch
- Manually adjustable backrest and armrests
- 3 heavy-duty castors allow for easy movement around the operation theatre
- Foot operated castor break and manual blocking system provide stability whenever required
- Battery pack - AC/DC power supply operated

Technical specification

Height adjustment range	550 - 750mm
Arm supports rotation range	up to 129°
Backrest height adjustment	75mm
Backrest forward/backward adjustment	55mm
Supply voltage*	~100 - 240V AC, 50/60Hz
Battery pack	24V, 2.2Ah
Battery charging	24V ±0.5V, 1.2A
Working time	up to 20h
Packaging	1 pallet; 120x80x160cm; 85kg (gross)

Dimensions

Width	885 - 1 015mm
Depth	740mm
Height	835 - 1 110mm



vanilla



mint



juniper



very berry

*- The power supply can support both 100VAC and 240VAC power sources. Voltage value detection build-in.