



AJL Ophthalmic, S.A. is the only Spanish company (July 1992) dedicated to the r&D, manufacture and marketing of medical devices in the field of Ophthalmology.

The company has reached a leading market postion thanks to the fulfillment of very exigent Regulatory and Quality Control Requirements, 3000 m2 facilities with cutting-edge tecnology as well as the most qualified personnel for the manufacture of its products.



AJL has an extensive sales network with presence in more than 80 countries.

AJL VET



 CANINE IOLS 10.25MM TXY65102S 11.25MM TXY65112S 12.25MM TXY65122S 13.25MM TXY65132S 14.25MM TXY65142S 15.25MM TXY65152S 41D 6.5 OPT <i>Disposable injector is supplied</i>	 FELINE IOLS 13MM CAT701300 53.5D 7.0 OPT <i>Disposable injector is supplied</i>	 PMMA CANINE IOLS 14MM - 6.5 OPT - P651400 16MM - 7.0 OPT - P701600 41D	 TRYPAN BLUE 0.6mg/ml BLUEVET 5x 1ML SYRINGES/BOX <i>Local regulations may apply.</i>
 CTR5 11MM ACV001102 13MM ACV001302 15MM ACV001502 13MM* ACV001302PL 16MM* ACV001602PL <i>*PRE-LOADED</i>	 VISCOELASTICS NaHA 1.4% VV14 NaHA 2% VV20 NaHA 3% VV30 HPMC 2% VCELL 5x 2ML SYRINGES/BOX	 SILICONE OIL 5000 cSt 55000VET 10ML SYRINGE	 MEMBRANES 2x2cm 100um MCV22100 2x2cm 200um MCV22200 <i>COLLAGEN TYPE I</i>



AJL Ophthalmic S.A.
Parque Tecnológico de Álava
C/ Ferdinand Zeppelin 1
01510 Miñano, Álava - Spain
T. +34 945 298 256
F. +34 945 298 209
ajlsa@ajlsa.com
www.ajlsa.com



The future of vision

AJL PRODUCT PORTFOLIO

AJL Ophthalmic S.A.
Parque Tecnológico de Álava
C/ Ferdinand Zeppelin 1
01510 Miñano, Álava - Spain
T. +34 945 298 256
F. +34 945 298 209
ajlsa@ajlsa.com
www.ajlsa.com

CORNEA



keratoconus




Intacs® SK (6mm optical zone)



AJL-RING

AFR (5 mm optical zone)
AFR6 (6 mm optical zone)



AJL pro+

CW (clockwise/5 mm optical zone)
CW6 (clockwise/6 mm optical zone)
CCW (counterclockwise/5 mm optical zone)
CCW6 (counterclockwise/6 mm optical zone)

VACUUM SYSTEM
VACUMM SYSTEM KV2000



BRIEFCASE INCLUDED
Dimensions: 46x33x15
Weight: 5,10 Kg.

Medical device designed to stabilize the cornea and facilitate the tunneling process during the intracorneal rings implantation surgery

Intercorneal Implant for Corneal Transplant



ENDO-K PRO

EQP75 7.50mm
EQP80 8.00mm
EQP85 8.50mm
EQP90 9.00mm



NEORING

AQ080 8.00mm
AQ085 8.50mm
AQ090 9.00mm



AJL-OOKP

Osteo-odonto-queratoprótesis

CATARACT



Iris expander



CANABRAVA
RING
CANA-RING

PMMA Intraocular lenses



AJL PMMA
P651300
(Posterior chamber)



AJL PMMA
A601250
(Anterior chamber)

Capsule stain solution



AJL BLUE
0,06%

Hydrophilic Intraocular lenses



AIALA
F601250



AS-IOL
LLASY60

Hydrophobic Intraocular lenses



AIALA DRY
LLASHP60



AIALA DRY
PRELODED
LLASHP60-PL



CAPSULAR TENSION RING
AC
ACS
AC-PL
ACS-PL

Open diameter	10mm	11mm	12mm	12mm	13mm	14mm
Close diameter	9mm	10mm	10mm	10,5mm	11mm	12mm

Viscoelastic solutions



AJL VISC
● AJL VISC 1%
● AJL VISC 1,4%
● AJL VISC 2%
● AJL VISC 3%
● AJL CELL 2%

GLAUCOMA



Scleral implant
ESNOPER
V2000



Uveoscleral implant
ESNOPER CLIP
VCLIP



NPDS KIT
NPDS surgical instruments
(non-penetrating deep sclerectomy)

RETINA



NPB
MACULAR BUCKLE

ILM stain solution



AJL BBG



SIOBAL
S1000
S2000
S5000

SIOBAL

OCULOPLASTY



OCUFIT SPHERICAL

6316-POI 14mm	6317-POI 20mm
6326-POI 16mm	8703-POI 21mm
6327-POI 18mm	6322-POI 22mm
8702-POI 19mm	8704-POI 23mm



OCULFIT CONICAL

80016-POI 16mm
80018-POI 18mm
80020-POI 20mm



OCULFIT SHEET

8438-POI ULTRA THIN
7210-POI THIN
6330-POI STANDARD



OCULAR SHIELD

PROCUL-S-AJL
PROCUL-M-AJL
PROCUL-L-AJL



SYMBLEPHARON RING

SF-S-AJL
SF-M-AJL
SF-L-AJL



PERFORATED
COMFORMER

CF-S-AJL
CF-M-AJL
CF-L-AJL



ALOS PALPEBRAL
IMPLANTS

GOLD TH (0,6 a 1,8)