



AS-IOL

Intraocular Lens

Brand name	AS-IOL
Reference	LLASY60
Description	One-piece posterior chamber intraocular lens
Optic design	Monofocal, neutral asphericity
Material	Hydrophilic acrylic
Filter	UV & Blue Light
Power Range	-5.0 D to +40.0 D in 1.0 D increments / +12.0 D to +25.0 D in 0.5 D increments
Dimensions	Optic diameter: 6.00 mm Overall diameter: 10.75 mm
Edge	360° Square
Haptic Design	Haptics in closed C-Loop
A-Constant / ACD	A-Constant: (optical) 118.0 A-Constant: (u/s) 117.3 ACD: 4.56
Water content	25.5 %
Recommended injection system	Viscoject 2.2
Sterilization	Steam
Shelf life	3 Years

CE 2797

ajlsa.com

