




AIALA DRY

Intraocular Lens

Brand name	AIALA DRY
Reference	LLASHP60-PL
Description	One-piece posterior chamber intraocular lens
Optic design	Monofocal, neutral asphericity
Material	Hydrophobic acrylic
Filter	UV & Blue Light
Power Range	-5.0 D to +40.0 D in 1.0 D increments / +10.0 D to +30.0 D in 0.5 D increments
Dimensions	Optic diameter: 6.00 mm Overall diameter: 13.00 mm
Edge	360° Square
Haptic Design	Haptics in open C-Loop
A-Constant / ACD	A-Constant: (optical) 119.7 A-Constant: (u/s) 119.2 ACD: 5.66
Water content	< 4 %
Injection System	Pre-loaded on Accuject 2.2 injector system
Sterilization	Ethylene Oxide (ETO)
Shelf life	3 years



CE 2797

ajlsa.com

