



AIALA

Intraocular Lens

Brand name	AIALA
Reference	F601250
Description	One-piece posterior chamber intraocular lens
Optic design	Monofocal spherical
Material	Hydrophilic acrylic
Filter	UV
Power Range	-10.0 D to +40.0 D in 1.0 D increments / +12.0 D to +25.0 D in 0.5 D increments
Dimensions	Optical diameter: 6.00 mm Overall diameter: 12.50 mm
Edge	360° Square
Haptic design	Haptics in open C-Loop
A-Constant / ACD	A-Constant: (optical) 118 A-Constant: (u/s) 117.30 ACD: 4.56
Water content	25.5 %
Recommended injection system	Viscoject 2.2
Sterilization	Steam
Shelf life	3 Years

CE 0318

ajlsa.com

