




AJL PMMA

Intraocular Lens

Brand name	AJL PMMA	
Reference	P651300 	A601250 
Description	One-piece posterior chamber intraocular lens	One-piece anterior chamber intraocular lens
Optic design	Spherical monofocal	Spherical monofocal
Material	Polymethylmethacrylate (PMMA)	Polymethylmethacrylate (PMMA)
Filter	UV	UV
Power range	-10.0 D to +40.0 D in 1.0 D increments / +12.0 D to +25.0 D in 0.5 D increments	+6.0 D to +30.0 D in 1.0 D increments
Dimensions	Optic diameter: 6.50 mm Overall diameter: 13.00 mm	Optic diameter: 6.00 mm Overall diameter: 12.50 mm
Edge	360° Square	360° Square
Haptic design	Haptics in open C-Loop	Haptics in open Z-Loop
A-Constant / ACD	A-Constant: (optical) 118.6 A-Constant: (u/s) 118.0 ACD: 4.96	A-Constant: (optical) 116.2 A-Constant: (u/s) 115.5 ACD: 3.51
Water content	0 %	0 %
Sterilization	Ethylene Oxide (ETO)	Ethylene Oxide (ETO)
Shelf life	4.5 years	4.5 years

CE 0318

